

Endoscopy Image Processing & Classification

Sujal, Al20BTECH11020

Artificial Intellegince

Indian Institute of Technology Hyderabad, Telangana, India

Table of Content

- Introduction
- Dataset
- Pre-processing
- Model Implementation
- Training
- Results
- Conclusion
- References



INTRODUCTION

- Diagnosing gastrointestinal diseases of images taken with the Wireless Capsule Endoscopy (WCE) method.
- Capture visual images of the GI tract.
- The classification of images taken using the WCE method, into four classes: polyps, ulcerative colitis, esophagitis, and normal
- Comparing EfficientNetB2 Model & InceptionNetv3 Model.

Dataset

Table 1. Dataset specifications.

| Raw samples | Class | Train | Validation | Test | Total |
|-------------|-------------|-------|------------|------|-------|
| | Normal | 800 | 500 | 200 | 1500 |
| | Ulcer | 800 | 500 | 200 | 1500 |
| | Polyps | 800 | 500 | 200 | 1500 |
| | Esophagitis | 800 | 500 | 200 | 1500 |
| Total | | 3,200 | 2,000 | 800 | 6,000 |

Dataset size: 1.2GB





Figure 2. Gastrointestinal conditions: polyp (a), ulcerative colitis (b), esophagitis (c), healthy gastrointestinal tract (d)



Pre-processing by using Keras' "ImageDataGenerator"

- Images are resized to 224 x 224 x 3 pixels.
- train(53.3%), validate(33.3%) and test(13.3%) images.

| Augmentation | Value | |
|-----------------|-------|--|
| Horizontal flip | True | |
| Rotation Range | 15 | |
| Shear Range | 0.2 | |
| Zoom Range | 0.2 | |
| Height shift | 0.1 | |
| Width shift | 0.1 | |





Model Implementation


- Using transfer learning by using pre-trained models on ImageNet.
- InceptionNetv3 and EfficientNetB2 to extract features from an image.
- Use the Adam optimizer due to it converge faster and more efficiently.
- Record all metrics like accuracy, Precision, Recall and AUC at each epoch.

InceptionNet Vs EfficientNet


| ayer (type) | Output Shape | Param # |
|--|--------------------|----------|
| inception_v3 (Functional) | (None, 5, 5, 2048) | 21802784 |
| flatten (Flatten) | (None, 51200) | 0 |
| dense (Dense) | (None, 512) | 26214912 |
| batch_normalization_94 (Bat chNormalization) | (None, 512) | 2048 |
| gaussian_noise (GaussianNoi se) | (None, 512) | 0 |
| dropout (Dropout) | (None, 512) | 0 |
| dense_1 (Dense) | (None, 4) | 2052 |
| Total params: 48,021,796 Trainable params: 26,217,988 Non-trainable params: 21,803 | , 808 | |

| Model: "sequential" | | | | |
|---|--------------------|-------------|--|--|
| Layer (type) | Output Shape | Param # | | |
| efficientnetb2 (Functional) | (None, 7, 7, 1408) | 7768569 | | |
| gaussian_noise (GaussianNoi se) | (None, 7, 7, 1408) | | | |
| global_average_pooling2d (G lobalAveragePooling2D) | (None, 1408) | | | |
| dense (Dense) | (None, 256) | 360704 | | |
| batch_normalization (BatchN ormalization) | (None, 256) | 1024 | | |
| gaussian_noise_1 (GaussianN oise) | (None, 256) | | | |
| dropout (Dropout) | (None, 256) | | | |
| dense_1 (Dense) | (None, 4) | 1028 | | |
| Total params: 8,131,325 | | | | |
| Trainable params: 362,244 | | | | |
| Non-trainable params: 7,769,081 | | | | |

Training:



InceptionNetV3 epochs: 40




EfficientNetB2 epoch: 20

Results using Test data

| | precision | recall | f1-score |
|--------------------------|-----------|--------|----------|
| normal cells | 0.70 | 0.97 | 0.81 |
| ulcerative colitis cells | 0.79 | 0.62 | 0.70 |
| polyps cells | 0.75 | 0.66 | 0.70 |
| esophagitis cells | 0.96 | 0.91 | 0.94 |
| Average | 0.80 | 0.79 | 0.79 |

Table 2. Classification Report of inceptionV3 model on Test data







Table 3. Classification Report of EfficientNetB2 model on Test data



Conclusion

- EfficientNetB2 with an accuracy of 94%, perform better than InceptionNetV3, with an accuracy of 79%.
- InceptionNetv3 uses a combination of convolutional layers, max-pooling layers, and fully connected layers.
- EfficientNetB2 uses a compound scaling method to balance the network's depth, width, and resolution.
- EfficientNetB2 requires less training time than InceptionNetv3 due to its efficient architecture and compound scaling method.
- EfficientNetB2 has a smaller model size compared to InceptionNetv3, making it easier to deploy on resource-constrained devices.

References

- KVASIR Dataset: https://dl.acm.org/doi/abs/10.1145/3083187.3083212
- Diagnosing gastrointestinal diseases from endoscopy images through a multifused CNN: https://www.sciencedirect.com/science/article/pii/S1746809422002051
- Endoscopic Image Classification Based on Explainable Deep Learning
 https://www.mdpi.com/1424-8220/23/6/3176
- Transfer Learning using Tensorflow: https://www.tensorflow.org/guide/keras/transfer_learning
- InceptionV3: https://www.tensorflow.org/api_docs/python/tf/keras/applications/inception_v3/Inception V3
- EfficientNetV2: https://www.tensorflow.org/api_docs/python/tf/keras/applications/efficientnet_v2/Efficient NetV2L

